

*Cleaned and Packaged in a Class 1000 Controlled Environment.

Meets ANSI/ESD S1.1 wristband resistance interior < 100 kilohms (<1 x 10^5 ohms) Meets ANSI/ESD S1.1 wristband resistance exterior >10 megohms (>1 x 10^7 ohms) Meets ANSI/ESD S20.20 and Wrist Strap standard ANSI/ESD S.1


Unless otherwise noted, tolerance ± 10% Specifications and procedures subject to change without notice.

- A. ABS cap with slotted opening to allow size adjustment.
- B. 300 Series Stainless Steel Backplate forms channel to hold up to 10 links of band allowing size adjustment.
- C. 4mm snap stud machined from 300 series stainless steel. Fits to 4mm spring loaded medical type socket. Coil cord disengages from snap with 1-5 lbs. of perpendicular force.
- D. Inside links constructed of 300 Series Stainless Steel. Width of links is 0.630".
- E. .002" polyphenylsulfone sheet is laminated to each link.Provides >5,000 volt RMS insulation per ASTM D149. Scratch resistant.
- F. Size is adjusted by sliding 1 to 10 links of the band into the backplate "C" channel.
- G. 300 Series Stainless Steel clip for electrical contact between snap stud and backplate.
- H. Underwriters Laboratory listing 90P1C.



United States of America

Premium Metal Expansion Wristband



DESCO EAST: ONE COLGATE WAY, CANTON, MA 02021-1407

PHONE (781) 821-8370 FAX (781) 575-0172

DRAWING NUMBER 09041 DATE: January 2016